



# Instructions For Use

## HCS1013-IFU

205 South 600 West Logan, Utah 84323, U.S.A. – Tel. (800) 729-8350 – Tel. (435) 755-9848 – Fax (435) 755-0015 – www.scytek.com Rev. 5, 5/16/2023

## Iron Control Slides

### Description

Iron Control Slides contain formalin fixed paraffin embedded sections of Human or Animal tissue cut at 4 micronthickness. Sections are known to contain deposits of ferric iron.

### Storage

Store slides at 2-25°C

### Suggested Stain Kit (not provided)

IRN-1

See below procedure

## Iron Stain Kit

### Description and Principle

The Iron Stain Kit is intended for use in the detection of ferric iron in tissues and bone marrow. Ferric iron is normally found in small amounts in bone marrow and the spleen. Abnormally large deposits may be seen in hemochromatosis and hemosiderosis.

This product is based on the Prussian Blue reaction in which ferric iron reacts with an acidic solution of potassium ferrocyanide to form an insoluble blue product called Prussian blue. Only ferric ions loosely bound to protein complexes will stain blue. Strongly bound ferric ions will not stain.

### Expected Results

Iron:	Bright Blue
Nuclei:	Red
Background:	Pink

### Smears

**Sideroblasts:** These are nucleated erythrocytes containing at least one small blue granule. If the blue granules surround the nucleus, the cell is a ringed sideroblast.

**Siderocytes:** These are non-nucleated erythrocytes containing at least one blue granule.

**Reticuloendothelial Iron:** Usually seen as blue particles on the marrow smear or as blue particles in the cytoplasm or phagocytic cells.

### Kit Contents

	<u>Storage</u>
1. Potassium Ferrocyanide Solution	18-25°C
2. Hydrochloric Acid Solution (2%)	18-25°C
3. Nuclear Fast Red Solution	18-25°C

### Suggested Controls (not provided)

Spleen, Bone Marrow.

### Uses/Limitations

For In-Vitro Diagnostic use only.

Do not use if reagents become cloudy or precipitate

Do not use past expiration date.


Use caution when handling reagents.

Non-Sterile

Intended for FFPE sections cut at 5-10µm.

This procedure has not been optimized for frozen sections.

Frozen sections may require protocol modification.



Ferric iron deposits in Human Liver viewed at 20X

### Storage

Store kit and all components at room temperature (18-25°C).

### Safety and Precautions

Please see current Safety Data Sheets (SDS) for this product and components GHS classification, pictograms, and full hazard/precautionary statements.

### Procedure

**Note:** Use acid-washed or bleach-washed glassware. Rinse all glassware with distilled water prior to use. Do not use metal forceps to transfer slide during staining procedure.

1. Deparaffinize sections if necessary and hydrate to distilled water.
2. Mix equal volumes of Potassium Ferrocyanide Solution and Hydrochloric Acid Solution to make a working Iron Stain Solution. Use once and discard.
3. Incubate slide in working Iron Stain Solution for 3-5 minutes.
4. Rinse thoroughly in distilled water.
5. Stain slide in Nuclear Fast Red Solution for 5 minutes.
6. Rinse in 4 changes of distilled water.
7. Dehydrate in 3 changes of absolute alcohol for 2 minutes each.
8. Clear, and mount in synthetic resin.

## References

1. Ryosuke Nakamura, Yasukazu Takanezawa, Yuka Ohshiro, Shimpei Uraguchi, Masako Kiyono, Effects of chemical forms of gadolinium on the spleen in mice after single intravenous administration, *Biochemistry and Biophysics Reports*, Volume 29, 2022, 101217, ISSN 2405-5808, <https://doi.org/10.1016/j.bbrep.2022.101217>.
2. Yamaguchi, H., Honda, S., Torii, S. et al. Wipi3 is essential for alternative autophagy and its loss causes neurodegeneration. *Nat Commun* 11, 5311 (2020). <https://doi.org/10.1038/s41467-020-18892-w>
3. Jing Qi, Jong-Won Kim, Zixiong Zhou, Chae-Woong Lim, Bumseok Kim, Ferroptosis Affects the Progression of Nonalcoholic Steatohepatitis via the Modulation of Lipid Peroxidation-Mediated Cell Death in Mice, *The American Journal of Pathology*, Volume 190, Issue 1, 2020, Pages 68-81, ISSN 0002-9440, <https://doi.org/10.1016/j.ajpath.2019.09.011>.
4. G.-H. Young, T.-M. Huang, C.-H. Wu, C.-F. Lai, C.-C. Hou, K.-Y. Peng, C.-J. Liang, S.-L. Lin, S.-C. Chang, P.-R. Tsai, K.-D. Wu, V.-C. Wu, and W.-J. Ko, "Hemojuvelin Modulates Iron Stress During Acute Kidney Injury: Improved by Furin Inhibitor," *Antioxidants & Redox Signaling*, vol. 20, no. 8, pp. 1181-1194, Mar. 2014.
6. Tetsuya Kagawa, Yuki Matsumi, Hiromichi Aono, Toshiaki Ohara, Hiroshi Tazawa, Kunitoshi Shigeyasu, Shuya Yano, Sho Takeda, Yasuhiro Komatsu, Robert M. Hoffman, Toshiyoshi Fujiwara & Hiroyuki Kishimoto (2021) Immuno-hyperthermia effected by antibody-conjugated nanoparticles selectively targets and eradicates individual cancer cells, *Cell Cycle*, 20:13, 1221-1230, DOI: 10.1080/15384101.2021.1915604
5. Sheenan, D.C., Hrapchak, B.B. *Theory and Practice of Histotechnology*, 2nd Edition. Battelle Press, Columbus, OH. Page 217. 1980
6. Carson, F.L., *Histotechnology: A Self-Instructional Text*, ASCP Press, Chicago, IL. Pages 214-215. 1990



ScyTek Laboratories, Inc.  
205 South 600 West  
Logan, UT 84321  
435-755-9848  
U.S.A.