

P.O. Box 3286 - Logan, Utah 84323, U.S.A. - Tel. (800) 729-8350 - Fax (435) 755-0015 - www.scytek.com

Amyloid Stain Kit (Congo Red)

Description: The Amyloid Stain Kit (Congo Red) is intended for use in the histological visualization of amyloid in tissue sections. Examination under a polarizing microscope results in green birefringence of amyloid.

Amyloid:	Red to Pink
Erythrocytes:	Light Orange
Eosinophil Granules:	Orange to Red
Nuclei:	Blue

Uses/Limitations: For In-Vitro Diagnostic use only.
 Histological applications.
 Do not use past expiration date.
 Use caution when handling these reagents.


Control Tissue: Freshly cut sections containing amyloid.
 Cut sections 6-12 microns to show smaller amyloid deposits.


Availability/Contents:

<u>Item #</u>	<u>Kit Contents</u>	<u>Storage Conditions</u>
CRA030	Congo Red Solution (30 ml) x2	Room Temperature
HMM030	Hematoxylin (30 ml)	Room Temperature
BRT030	Bluing Reagent 30 ml)	Room Temperature

Required but not included: 95% Ethyl Alcohol
 100% Ethyl Alcohol

Precautions: Congo Red Solution is flammable.
 Keep away from open flame.
 Avoid contact with skin and eyes.
 Harmful if swallowed.
 Follow all Federal, State, and local regulations regarding disposal.
 Use in chemical fume hood whenever possible.

Storage: 18° C  25° C

 ScyTek Laboratories, Inc.
 205 South 600 West
 Logan, UT 84321
 435-755-9848
 U.S.A.

CE Authorized Representative in Europe
 (Regulatory affairs only)
 EmergoEurope (31)(0) 70 345-8570
 Molsnstraat 15
 2513 BH Hague, The Netherlands

Instructions For Use

AMY-2-IFU

Rev. Date: August 2, 2017

Revision: 3

Page 2 of 2


P.O. Box 3286 - Logan, Utah 84323, U.S.A. - Tel. (800) 729-8350 - Fax (435) 755-0015 - www.scytek.com


Procedure (Standard):

1. Deparaffinize sections if necessary and hydrate to distilled water.
2. Apply 4-6 drops of Hematoxylin and incubate for 5 minutes.
3. Rinse slide in tap water.
4. Apply 4-6 drops of Bluing Reagent and incubate for 30 seconds.
5. Rinse slide in distilled water.
6. Dip slide in 95% alcohol for 5 seconds.
7. Pour one bottle of Congo Red Solution into provided coplin jar and incubate slide for 20 minutes. **Note:** Ensure that tissue is fully immersed in the coplin jar.
8. Dip twice (quickly) in 100% alcohol.
9. Dip repeatedly (4-5 dips) in clearing agent, and mount in synthetic resin.

References:

1. Puchtler, H, et al: On the binding of Congo Red amyloid. J. Histochem. Cytochem. Vol. 10: pages 355-363, 1962.
2. Eastwood, H. & Cole, K.R., Staining of amyloid in buffered Congo Red in 50% ethanol. Stain Technology. Vol. 46: pages 208-209, 1971.
3. Carson, F.L., Histotechnology; A Self-Instructional Text, 2nd Edition. ASCP Press, Chicago, IL. Pages 117-121, 1996.
4. Churukian, C., Improved Puchtler's Congo Red method. J. of Histotechnology. Vol. 23: pages 139-141, 2000.

Storage: 18° C  25° C



ScyTek Laboratories, Inc.
205 South 600 West
Logan, UT 84321
435-755-9848
U.S.A.

 Authorized Representative in Europe
 (Regulatory affairs only)
EmergoEurope (31)(0) 70 345-8570
Molsnstraat 15
2513 BH Hague, The Netherlands