



Instructions For Use

BRT-IFU

205 South 600 West Logan, Utah 84323, U.S.A. – Tel. (800) 729-8350 – Tel. (435) 755-9848 – Fax (435) 755-0015 – www.scytek.com Rev. 5, 5/21/2024

Bluing Reagent for Hematoxylin

Description and Principle

A pH-controlled solution for quick and effective bluing of hematoxylin used in histological and cytological procedures.

Expected Results

Nuclei: Purple
Nuclei after Bluing: Blue

Additional Required Reagents (Not Included)

1. Hematoxylin, Mayer's (Lillie's Modification)

Suggested Controls (not provided)

Any well-fixed tissue stained with hematoxylin.

Uses/Limitations

For In-Vitro Diagnostic use only.
Do not use if reagents become cloudy or precipitate
Do not use past expiration date.
Use caution when handling reagents.
Non-Sterile
Intended for FFPE sections cut at 5-10µm.
This procedure has not been optimized for frozen sections.
Frozen sections may require protocol modification.

Storage

Store at room temperature (18-25°C).

Safety and Precautions

Please see current Safety Data Sheets (SDS) for this product and components GHS classification, pictograms, and full hazard/precautionary statements.

Procedure

Counterstaining for IHC:

Dip slide in Mayer's Hematoxylin (Lillie's Modification) several times. Blue in bluing reagent for 15-30 seconds.

-OR-

Incubate in Mayer's Hematoxylin (Lillie's Modification) for 30-60 seconds. Blue in bluing reagent for 15-30 seconds.

H&E staining and Hematoxylin standalone:

1. Stain for 3-5 minutes in Mayer's Hematoxylin (Lillie's Modification). Note: Longer incubation times provide a darker stain.
2. Blue in bluing reagent for 15-30 seconds.
3. Continue with incubation in Eosin or dehydration and clearing and mounting.




Fig 1. IHC Staining with CD34 on Human tonsil. Counterstained with Mayer's Hematoxylin (Lillie's Modification) for 30 seconds followed by bluing with Bluing Reagent.





Fig 2. Blued Hematoxylin on Pig Lung.

References

1. Hanwright PJ, Rath JB, von Guionneau N, Slavin B, Pinni S, Zlotolow D, Shores J, Dellon AL, Tuffaha SH. The effects of a porcine extracellular matrix nerve wrap as an adjunct to primary epineurial repair. *The Journal of Hand Surgery*. 2021 Sep 1;46(9):813-e1.
2. Jacobs SO, Sheller-Miller S, Richardson LS, Urrabaz-Garza R, Radnaa E, Menon R. Characterizing the immune cell population in the human fetal membrane. *American Journal of Reproductive Immunology*. 2021 May;85(5):e13368.
3. Albracht-Schulte K, Gonzalez S, Jackson A, Wilson S, Ramalingam L, Kalupahana NS, Moustaid-Moussa N. Eicosapentaenoic Acid Improves Hepatic Metabolism and Reduces Inflammation Independent of Obesity in High-Fat-Fed Mice and in HepG2 Cells. *Nutrients*. 2019 Mar;11(3):599.
4. Bellon M, Moles R, Chaib-Mezrag H, Pancewicz J, Nicot C. JAG1 overexpression contributes to Notch1 signaling and the migration of HTLV-1-transformed ATL cells. *Journal of Hematology & Oncology*. 2018 Dec;11(1):1-2.
5. Nandi SS, Shahshahan HR, Shang Q, Kutty S, Boska M, Mishra PK. MiR-133a mimic alleviates T1DM-induced systolic dysfunction in Akita: an MRI-based study. *Frontiers in physiology*. 2018 Oct 10;9:1275.
6. Jung SH, Song HY, Hyun YS, Kim YC, Whang I, Choi TY, Jo S. A brain atlas of the long arm octopus, *octopus minor*. *Experimental neurobiology*. 2018 Aug;27(4):257.

7. Bellon M, Lu L, Nicot C. Constitutive activation of Pim1 kinase is a therapeutic target for adult T-cell leukemia. Blood, The Journal of the American Society of Hematology. 2016 May 19;127(20):2439-50.

 ScyTek Laboratories, Inc.
205 South 600 West
Logan, UT 84321
U.S.A.

Emergo Europe
Westervoortsedijk 60
6827 AT Arnhem, The Netherlands