



# Instructions For Use

## SMG-IFU

205 South 600 West Logan, Utah 84323, U.S.A. – Tel. (800) 729-8350 – Tel. (435) 755-9848 – Fax (435) 755-0015 – www.scytek.com Rev. 3, 5/23/2024

### Mucicarmin Solution (Southgate's)

#### Description and Principle

Mucicarmin Solution (Southgate's) is a component of the Mucicarmin Stain Kit (Catalog# SMS-1) and is intended for use in the histological visualization of acid mucopolysaccharides in tissue sections. Mucicarmin Solution is the component responsible for staining Mucin red. This product is useful in distinguishing mucin negative undifferentiated squamous cell lesions from mucin positive adenocarcinomas. In addition, this product will stain the mucopolysaccharide capsule of *Cryptococcus neoformans*.

#### Expected Results

Mucin:	Pink/Red
Capsule of <i>Cryptococcus</i> :	Red
Nuclei:	Blue to Green
Other tissue components:	Yellow

#### Kit Contents

#### Storage

#### Additional Kit Reagents Sold Separately

1. Hematoxylin, Mayer's	18-25°C
2. Bluing Reagent	18-25°C
3. Mucicarmin Solution	2-8°C
4. Tartrazine Solution	18-25°C

#### Suggested Controls (not provided)

Colon, Intestine, Bronchi Epithelial cells

#### Uses/Limitations

For In-Vitro Diagnostic use only.

Do not use if reagents become cloudy or precipitate

Do not use past expiration date.

Use caution when handling reagents.

Non-Sterile

Intended for FFPE sections cut at 5-10µm.

This procedure has not been optimized for frozen sections.

Frozen sections may require protocol modification.

#### Storage

Store at 2-8°C.

#### Safety and Precautions

Please see current Safety Data Sheets (SDS) for this product and components GHS classification, pictograms, and full hazard/precautionary statements.

#### Procedure

1. Deparaffinize sections if necessary and hydrate to distilled water.
2. Stain slide in Hematoxylin for 2-3 minutes.
3. Rinse for 2 minutes in running tap water.
4. Dip slide in Bluing Reagent for 30 seconds.
5. Rinse slide in distilled water.
6. Stain slide in Mucicarmin Solution for 10 minutes.
7. Rinse quickly in tap water followed by a quick rinse in distilled water.
8. Stain slide in Tartrazine Solution for 1 minute.




Epithelial mucins in Human GI tract demonstrated with Mucicarmin Stain Kit

9. Rinse slide 1 time in absolute alcohol.
10. Dehydrate slide quickly in 3 changes of absolute alcohol.
11. Clear in 2 changes of xylene or xylene substitute, and mount in synthetic resin.

#### References

1. Sheenan, D.C., Hrapchak, B.B. Theory and Practice of Histotechnology, 2nd Edition. Battelle Press, Columbus, OH. Pages 172-173.

 ScyTek Laboratories, Inc.  
205 South 600 West  
Logan, UT 84321  
U.S.A.

 Emergo Europe  
Westervoortsedijk 60  
6827 AT Arnhem, The Netherlands