

Polyclonal Antibody to Human Alpha-1-Fetoprotein (AFP)

Description:

- Immunogen:** Alpha-1-fetoprotein isolated from human cord serum.
- Format:** This antibody has been pretitered and quality controlled to work on formalin-fixed paraffin-embedded as well as acetone fixed cryostat tissue sections. No further titration is required.
- Specificity:** This antibody reacts with alpha-1-fetoprotein. A positive staining in the hepatocytes of fetal tissue was observed. This antibody does not react with any other cell or tissue types.

Uses/Limitations: Immunohistochemistry
For Research Use Only.
Do not use past expiration date.

Storage: 2-8° Centigrade.

Procedure: We suggest an incubation period of 30-60 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation should be determined by the user.

Precautions: Contains Sodium Azide as a preservative.
Do not pipette by mouth.

References: i) MacDonald et al. Clin Chim Acta 87: 367, 1978.
ii) Jacobson et al. Am J Surg Pathol 5: 257, 1981.